

CONTINUOUS, REAL-TIME CARDIOVASCULAR MONITORING



Beat to Beat

Stroke Volume Variation

LiDCO

Cardiac Sensor Systems

Introducing the LiDCO^{plus} Hemodynamic Monitor

CONTINUOUS, REAL-TIME MONITORING OF CARDIAC OUTPUT


Monitoring of the key cardiovascular parameters of blood pressure, cardiac output and oxygen delivery is essential for many major surgical and acutely ill patients. In addition, there are many other patients that could benefit from real-time cardiovascular monitoring if it were available in a safe and easy to apply manner. This was the philosophy that resulted in the development, design and manufacture of the unique LiDCO^{plus} Hemodynamic Monitor.

The LiDCO^{plus} system is a combination of the two innovative and novel monitors - the LiDCO System indicator dilution cardiac output monitor and the PulseCO System real time, continuous arterial waveform monitor, produced by LiDCO Ltd.

This unique combination provides beat-to-beat measurement of cardiac output with lower risk and high precision ⁶⁸.

INTRODUCING LiDCO LTD

- LiDCO Ltd's starting point for its products is with the requirements of the clinician, nurse and hospital patient for advanced minimally invasive cardiovascular monitoring.
- LiDCO Ltd's products result from the consolidation of a physiological knowledge base and leading physiological software and sensor technology.
- LiDCO Ltd's minimally invasive monitoring products improve hospital standards of care, improve patient outcome and reduce costs ⁶⁸.
- LiDCO Ltd's products provide medical staff with essential hemodynamic information, not just data, in an 'easy to interpret' manner.




The LiDCO^{plus} System combines the LiDCO & PulseCO Systems software and provides a real-time and continuous assessment of a patient's hemodynamic status.

The LiDCO™ System cardiac output method provides a bolus indicator dilution method of measuring cardiac output. A small dose of lithium chloride is injected via a central or peripheral venous line; the resulting arterial lithium concentration-time curve is recorded by withdrawing blood past a lithium sensor attached to the patient's existing arterial line. In terms of accuracy, clinical studies have demonstrated that over a wide range of cardiac outputs the LiDCO method is at least as accurate as thermodilution and even in patients with varying cardiac outputs⁷⁻¹². In one study LiDCO and thermodilution cardiac output were compared with an electromagnetic flow probe. The results of this study indicated that LiDCO had a higher precision compared with conventional bolus thermodilution cardiac output¹². The dose of lithium needed (0.15 - 0.3 mmol for an average adult) is very small and has no known pharmacological effects⁷⁴.

The PulseCO™ System software (incorporated in the LiDCO™*plus* Monitor) calculates continuous beat-to-beat cardiac output by analysis of the arterial blood pressure trace following calibration with an absolute LiDCO cardiac output value. This method has been shown to be accurate and reliable in various clinical settings. It has also been shown that recalibration is unnecessary for at least eight hours^{2-6, 79, 80} and more recently for 24 hours⁹⁴.

PULSECO SYSTEM AUTOCORRELATION ALGORITHM




The analogue arterial blood pressure trace is slaved from the conventional blood pressure monitor and undergoes a three step transformation


- Step 1 Arterial pressure transformation into a volume-time waveform (incorporating arterial tree compliance)
- Step 2 Deriving nominal stroke volume and heartbeat duration
- Step 3 Actual stroke volume via calibration with an absolute cardiac output value

Components of the LiDCO^{plus} Monitor

The PulseCO software calculates continuous beat-to-beat cardiac output by analysis of the arterial blood pressure trace following calibration with an absolute cardiac output value. This absolute cardiac output value is accurately and precisely measured using the innovative LiDCO lithium chloride bolus indicator dilution method.




The PulseCO software calculates the pulse power and derived stroke volume from the arterial waveform. This avoids the necessity for detection of any particular waveform features such as the dichrotic notch. Furthermore, arterial wave reflection does not have to be estimated. Due to this innovative and patented method of calculation, the PulseCO remains accurate and reliable over a wide range of hemodynamic states in surgical, post operative and intensive care settings². Studies have demonstrated that re-calibration is unnecessary for at least eight hours²⁻⁶ and more recently 24 hours⁹⁴.



The LiDCO^{plus} serves as a reliable alternative to continuous cardiac output monitoring with the pulmonary artery catheter and can easily be used intraoperatively, in the ICU/HDU, trauma or burns unit, and cath lab.

LiDCO^{plus}'s accuracy is ensured with the proven LiDCO^{plus} lithium indicator dilution calibration procedure that uses existing venous and arterial access, making it fast, cost effective and minimally invasive.



The LiDCO^{plus} calculates cardiac output continuously by analysis of the arterial blood pressure trace following calibration with the absolute LiDCO cardiac output value. The concept of estimating cardiac output from the arterial pressure waveform has been extensively researched with the first researchers (Erlanger and Hooker) publishing in 1904¹.

The LiDCO System is an innovative bolus indicator dilution method of measuring cardiac output and Intra-Thoracic Blood Volume (ITBV) - the innovation is the use of lithium chloride as the indicator. A small dose of lithium chloride is injected via a central or peripheral venous line (Fig 1); the resulting arterial lithium concentration-time curve is recorded by withdrawing blood past a lithium sensor attached to the patient's existing arterial line (Fig 2); the Monitor then calculates the cardiac output from the area of the primary dilution curve (Fig 3). The mean transit time of the lithium is derived for calculation of the ITBV.

The signal to noise ratio and hence accuracy for lithium is better than that seen with thermodilution - due to the fact that the lithium dose can be scaled to the size and cardiac output of the patient. Thermal noise from fluid infusion, respiration and patient warming has little, if any, effect on the lithium curve. The precision of the LiDCO System method means that only one lithium injection is required to accurately determine the cardiac output. In terms of accuracy, clinical studies have demonstrated that the single bolus LiDCO System method is at least as accurate as triplicate bolus thermodilution over a wide range of cardiac outputs and even in patients with varying cardiac outputs ⁷⁻¹².

The lithium chloride indicator dilution method of measuring cardiac output and intra-thoracic blood volume (ITBV)

Fig 2





Fig 1




Fig 3




Arterial Pressure

The arterial pressure waveform is slaved from your conventional blood pressure monitor via an analogue cable link. At the touch of a button the arterial pressure waveform is displayed on the LiDCO_{plus} screen. This can be viewed on the screen whenever the Trend, Graph and Chart screens are in use.

This Blood Pressure Waveform window is used to check the patient's systolic, diastolic and mean pressure. The pressure waveform shape and values should equate to those displayed on the primary blood pressure monitor.



This window also provides you with access to preload response values or volume status indicators of: Systolic Pressure Variation (SPV), Pulse Pressure Variation (PPV%), Stroke Volume Variation (SVV%) as well as Heart Rate Variation (HRV%).



$$SVV \% = (SV \text{ max} - SV \text{ min}) / [(SV \text{ max} + SV \text{ min}) / 2] \times 100$$


Patient with SVV less than 10% are unlikely to be preload responsive.

DYNAMIC VOLUME STATUS INDICATORS FOR VENTILATED PATIENTS

Thoracic pressure changes caused by mechanical ventilation induces cyclic changes in Left Ventricular Stroke Volume. These changes can provide an indication of the patient's ventricular preload status (see figure of Frank-Starling curve)^{18, 20, 26, 66, 81-83}.

By superimposing the Preload Response Parameters Window onto the Trend, Graph, or Chart screens, a continuous measurement of Stroke Volume, Systolic Pressure and Pulse Pressure Variation is displayed numerically and graphically.

For closed chest ventilated patients these volume status measurements provide a way of predicting fluid volume status and likely response to volume infusions. A fluid imbalance can have an adverse effect on a patient's cardiac performance and, in turn, oxygen delivery to key organs^{13-35, 50, 56, 66, 67, 81}.




Frank-Starling Curve: Patient Preload Status

Event Response

PRE-LOAD RESPONSIVENESS VIA THE EVENT RESPONSE DISPLAY

The Event Response display allows the user to view up to 2 hemodynamic variables in a higher resolution during a specific period (e.g. fluid challenge, inotrope change). The LiDCO_{plus} will also display percent change from start for each variable as a numeric value. Trend lines can also be added to the graphical display.

This feature is very useful when evaluating the patient's response to targeted interventions such as a fluid challenges or changes in inotrope therapy. The advantage this screen brings is in averaging the data display to smooth any noise and magnifying the display of changes so that they are easily viewed at a glance from a distance.



Benefits to you and your patient

BENEFITS OF LIDCO_{plus} CONTINUOUS CARDIAC OUTPUT MONITORING INCLUDE:

- Provides early warning of patient deterioration
- Optimisation of oxygen delivery
- Optimisation of fluid management
- Rational drug administration (e.g. Inotropes)
- The patient's condition is clearly communicated to clinical staff
- Reduces the work of health care staff
- Decrease the procedural complications
- Is minimally invasive and therefore widely applicable
- Is accurate
- Can be nurse driven
- Provides real time, beat-to-beat cardiac output and oxygen delivery
- Provides real time preload and afterload values
- Provides indexed values
- Provides easy data interpretation
- Provides bedside information management
- Has easy to use event markers
- Provides information not simply data
- Can be linked to most commonly found BP monitors
- Records historical data

BENEFITS OF THE LIDCO – LITHIUM CHLORIDE INDICATOR DILUTION CARDIAC OUTPUT METHOD:


- Provides an absolute cardiac output value via a novel and proven indicator dilution technique
- Provides ITBV
- Requires no additional invasive catheters to insert into the patient
- Is safe – using non-toxic bolus dosages
- Is simple and quick to set up
- Can be used with a range of LiCl dosages
- Is as accurate as meaned triplicate thermodilution
- Is not temperature dependent
- Is less invasive monitoring
- Utilises existing peripheral or central venous and arterial lines
- Can be set up and used by nursing staff
- Is a well studied and validated technique

The Trend Screen

A CLEAR CONTINUOUS DISPLAY


The Trend screen provides a continuous record of the patients hemodynamics. This large, clear and continuously updating screen displays the actual and indexed values for Cardiac Output, Systemic Vascular Resistance, Mean Arterial Pressure, Heart Rate and Left Ventricle Stroke Volume. Oxygen Delivery and Venous Oxygen Saturation can also be displayed if selected. The clinician can accurately track the patients' trend over several hours or minutes.

This screen facilitates intraoperative patient management allowing assessment of the immediate response to fluid challenge, drugs or other therapeutic interventions.



EVENT MARKER


Mark any number of events, such as initiation of inotropic agent or start of fluid infusion, using the event marker button. There are several standard descriptions that can be selected and additional detail added.



The marked event will appear as a small flag on the Trend and History screens and will be recorded in the LiDCOplus data file as the patient data is downloaded from the history screen. The event history window will allow selection of any events from the past 24 hours and take the user to the area on the hemodynamic history screen where the event was marked.

The Chart Screen

PARAMETER RELATIONSHIPS AT A GLANCE



$$\text{Pressure} = \text{Flow} \times \text{Resistance}$$

Showing the relationship between pressure, flow and resistance in an integral bar chart display. This screen simplifies the recognition and diagnosis of hemodynamic imbalance at the bedside. This means that potentially complex hemodynamic data can be easily interpreted and the necessary corrective actions taken quickly.


OXYGEN DELIVERY TARGETS AND $\text{ScvO}_2/\text{SvO}_2$ DISPLAY

Oxygen delivery optimisation has been demonstrated to be a key factor in the reduction in both length of stay and complications for post surgical patients in the ICU⁸⁸. Oxygen delivery targets can now be set in the patient limits. These will be used on both the Graph and Chart screens to enhance the display of data. Venous Oxygen Saturation can now be inputted to allow tracking as well as calculation of estimated oxygen consumption.

OXYGEN DELIVERY

The goal of monitoring cardiac output is to maximise the delivery of oxygen to the tissue beds. The LiDCO $plus$ monitor displays the parameter of oxygen delivery (and oxygen delivery index) in real time. The ability to have these two parameters monitored simultaneously on a real-time basis by the LiDCO $plus$ can have a major impact on patient care and outcomes. There is mounting evidence that monitoring oxygen delivery and cardiac output in at-risk patients can significantly reduce mortality and length of hospital stay⁵⁰.

In a recent study utilizing the LiDCO $plus$ it was shown that optimising oxygen delivery index (DO $_2$ I) to a target of 600ml/min/m 2 reduced morbidity by 50% and mean length of stay per patient by 12.3 days⁸⁸. Use of the LiDCO $plus$ facilitates the perioperative optimisation of patients.



HEMODYNAMIC PRE-OPTIMISATION IN HIGH-RISK PATIENTS


'High risk' surgery patients subjected to a reduction in global oxygen delivery are known to have increased levels of morbidity and mortality. Increasing global oxygen delivery has been reported to result in a dramatic improvement in outcome in these patients³⁶⁻⁶⁷. For example in the USA there are over 30 million operations performed annually, 10-15% of which (approximately 3 million operations) are deemed to be 'high risk'. These 'high risk' patients have an increased risk of death. A recent review of 21 randomised controlled trials with various approaches to treatment revealed statistically significant mortality reductions when patients with acute critical illness were treated early to achieve optimal goals before the development of organ failure⁵⁸. There is convincing evidence that measurement and manipulation of cardiac output, and therefore oxygen delivery, in selected patients reduces the risk of mortality^{59, 68, 69, 70}.

MORTALITY OUTCOMES FOLLOWING GOAL DIRECTED STUDIES:

Type of Surgery	Study Author (Year)	Mortality of Control Group (%)	Mortality of Goal Directed Treatment Group (%)
Vascular	Shultz et al (1985) ⁴⁴	29.0	2.9
General	Shoemaker et al (1988) ⁴⁵	33.0	4.0
Vascular	Berlanki et al (1991) ⁴⁶	9.5	1.5
Trauma	Fleming et al (1992) ⁴⁷	44.0	24.0
General & Vascular	Boyd et al (1993) ⁶⁷	22.2	5.7
Trauma	Bishop et al (1995) ⁴⁸	37.0	18.0
Hip Fracture	Sinclair et al (1997) ⁵¹	10.0	5.0
Peripheral Vascular	Ziegler et al (1997) ⁸⁴	9.0	5.0
Elective General	Wilson et al (1999) ⁵⁰	17.0	3.0
Elective Cardiac	Polonen et al (2000) ⁸⁵	3.0	1.0
General & Vascular	Lobo et al (2000) ⁸⁶	50.0	15.7

The Graph Screen

KEEPING PATIENT CARE ON TARGET




Allows easy assessment of the patients status at a glance and from a distance from the Monitor - useful in busy ICU's or high dependency units, as an early warning system to avoid adverse events. The graph screen is also ideal for the implementation of goal directed protocols that target oxygen delivery. The graph screen is the ideal "bedside" mode providing clear, easy-to-read feedback on changes to pressure, flow, and resistance. A patient specific "target zone" for Cardiac Output or Oxygen Delivery, Mean Arterial Pressure and thereby Systemic Vascular Resistance is defined specifically for each patient. The continuous display of the patient's last 12 heartbeats provides a ready reference to how well the therapeutic/hemodynamic targets are being maintained. The LiDCO_{plus} Monitor provides point of care information, not just more patient data.

CONTINUOUS, REAL-TIME MONITORING OF CARDIAC OUTPUT

The History Screen

KEEPING PATIENT RECORDS, FOR AUDIT, TEACHING, ANALYSIS AND RESEARCH



The History Screen can be used to look back over the last 24 hours of the patient's hemodynamic data. The display is similar to the display of the Trend Screen. Touch anywhere on the screen to see the parameter readings, at that time, or track back through the hours of data as required.

This screen is designed to aid in data collection: for clinical studies or simply to have a complete record of the patient's treatment you can record the patient's critical parameters at the touch of a button. The beat-to-beat patient data, event markers and calibrations are recorded either as a file for LiDCOview or as an Excel file on a USB memory stick. This powerful tool for data collection provides you with the ability to review, research and train using Historical data. Record several patients or several days of a single patient onto one easy to use USB memory stick.

A picture can be taken of the entire history screen. This can be stored on the hard drive or immediately downloaded to a USB memory stick as a JPEG file.

PulseCO Validation

A number of studies have been completed in various centres in Europe and the USA demonstrating the accuracy and precision of the PulseCO software (hosted within the LiDCO*plus* Hemodynamic Monitor) when used in cardiothoracic surgery, major surgery and the ICU ^{2-6, 79, 80, 94}.

INTRA OPERATIVE PERFORMANCE

One study at the University of Texas evaluated the change/drift in calibration factor across a post operative period of 8 hours in a group of 20 cardiac surgery patients ⁸⁰.

The results of the study found:


- that the range of cardiac outputs was 3.33 to 8.47 litres per minute;
- that, once calibrated, the PulseCO tracked the cardiac output continuously through the post surgery period without changes in calibration factor;
- that no significant differences were noted between the LiDCO, the thermodilution control and the PulseCO reading.

Differences in CO measurements by technique (Liters/min).

	0hours	2hours	4hours	6hours	8hours	all times
TLdiff	-0.1±0.1	0.0±0.1	0.1±0.1	0.0±0.2	-0.2±0.2	0.0±0.1
PLdiff	0.0±0.0	0.0±0.0	0.0±0.0	0.0±0.0	0.0±0.0	0.0±0.0
PTdiff	0.1±0.1	0.0±0.1	0.0±0.1	0.0±0.2	0.2±0.2	0.0±0.1

Conclusion:

'This technique appears to offer a safe, reliable and less invasive alternative to the traditional PA catheter for cardiac output monitoring in the immediate post-operative period after surgery.'





A study in intensive care patients ⁹⁴ has shown good correlation of LiDCO and PulseCO without recalibration over 24hrs.

POST OPERATIVE PERFORMANCE AND THE EFFECT OF DAMPING

A further study was undertaken at Duke University to examine the effects of arterial damping in a group of 22 patients in a surgical intensive care unit². The study was designed to evaluate the change/drift in calibration factor across a 24 hour period and in addition to examine whether the dynamic response (damping) of the arterial pressure monitoring system has an effect on the accuracy of the PulseCO software.

The results of the study found:

- that the range of cardiac outputs was 3.45 to 10.47 litres per minute;
- that there were acceptable bias and limits of agreement throughout the study with a correlation of $r^2 = 0.89$;
- that 68% of catheters had an undesirable dynamic response but that this did not affect the agreement between the measurements of cardiac output obtained by the PulseCO software and the Lithium dilution ($p = 0.976$)



Conclusion:

'Technical limitations of arterial catheter monitoring systems, such as low natural frequency or under or over damping, do not appear to influence the accuracy of the PulseCO measurement.'

'The pulse contour algorithm appears accurate despite a wide variety of arterial pressure waveform contours that are seen in clinical practice.'

'Re-calibration of the PulseCO software is recommended every four hours, but we have extended the period between calibrations to eight hours and demonstrated that the PulseCO software is still accurate after this time.'


The LiDCO^{plus} - Lithium Indicator Dilution

The LiDCO^{plus} bolus lithium indicator dilution method is used to provide an accurate actual cardiac output value for the patient. The lithium indicator dilution method is a very accurate and minimally invasive indicator dilution cardiac output measurement method. This is used to calibrate the PulseCO arterial waveform stroke volume value. It can also be used to calculate intra thoracic blood volume (ITBV).

THE PRINCIPLE

The bolus indicator dilution method of measuring cardiac output was first described by Henriques and developed by Hamilton et al in 1932⁷¹. This became widely adopted - the original technique using indocyanine green (ICG) as the marker. However, as this technique required frequent blood sampling and manual analysis of the dilution curve, it proved to be technically difficult and time consuming. The use of Lithium as an indicator of cardiac output was first described in 1993⁷² and has since been extensively validated⁷³. Lithium like ICG, is a non-diffusible indicator⁹². The method of using bolus indicator dilution to measure volume was described by Stewart⁹³ and the method of ITBV calculation is simply $ITBV = CO \times MTt$. MTt is the mean transit time of the lithium indicator from injection to detection.

The LiDCO^{plus} Monitor provides a lithium bolus indicator dilution method of measuring cardiac output. A small dose of lithium chloride is injected via a central or peripheral venous line; the resulting arterial lithium concentration-time curve is recorded by withdrawing blood past a lithium sensor.



ADVANTAGES OF THE LiDCO_{plus} LITHIUM INDICATOR DILUTION METHOD

The advantages of the LiDCO_{plus} method are that it is safe, accurate and simple to use:

- **Safe** - Central/peripheral venous and arterial catheters are usually already in place in patients needing cardiac output measurements. No further catheter is needed, so the method avoids the risks associated with pulmonary artery catheterisation. The method requires withdrawal of approximately 5 ml blood per determination; for an adult this is an insignificant amount. The injectate is a solution of lithium chloride. The dose needed (0.15 - 0.30 mmol for an average adult) is very small and has no known pharmacological effect ⁷⁴. The dosage regimen recommended is very conservative, making worst case assumptions on volume of distribution of lithium, patient weight (assumes 40kg) and absence of renal function.
- **Accurate** - Clinical trials have been completed that demonstrate that the LiDCO System is at least as accurate as thermodilution ^{11,12}.
- **Simple to use** - The method is simple and quick to use. It has the advantage that there is no unpleasant procedure for a conscious patient to undergo (such as insertion of a pulmonary artery catheter) and the time taken to set up and apply is between 5 and 10 minutes ⁶⁸.

THE SENSOR


The sensor consists of a lithium-selective electrode in a flow-through cell. It is disposable, sterilised by gamma irradiation and foil packed. The sensor is connected to a three-way tap on the arterial line and a small peristaltic pump restricts the flow through it to 4.5 ml/min. The flow-through cell is made of polycarbonate and designed with an eccentric inlet so that the blood swirls past the tip of the electrode.



The lithium-selective electrode is made of polyurethane with a central lumen. Silver chloride paint coats the inside and the outside. A wick, which is soaked in saline when the cell is first primed, makes the electrical connection between the blood in the cell in the vicinity of the tip of the electrode and the remote reference. This arrangement ensures adequate constancy of the voltage of the reference which is far enough from the blood to avoid a significant temperature effect. The electrode is filled with a reference material which provides a constant ionic environment and supports the membrane which is dip cast. The membrane is made of polyvinyl chloride and contains a lithium ionophore to make it selectively permeable to lithium ions.

Indicator Dilution Screen

GREATER PRECISION THAN A SINGLE BOLUS THERMODILUTION



DERIVATION OF ABSOLUTE CARDIAC OUTPUT VALUE

The voltage across the sensor membrane is related via the Nernst equation to the plasma $[Li^+]$. A correction is applied for plasma sodium concentration because in the absence of lithium the baseline voltage is determined by the sodium concentration. The voltage is measured using an amplifier optically isolated from the patient, then digitised on-line and analysed by the LiDCO^{plus} Monitor software.

Indicator dilution curves recorded in arterial blood consist of primary and secondary curves due to the initial circulation and then re-circulation of the indicator. Cardiac output is calculated as:

$$\text{Cardiac Output} = (\text{Lithium Dose} \times 60) / (\text{Area} \times (1 - \text{PCV}))$$

Where lithium dose is in mmol; Area is the integral of the primary curve (mM.s); PCV is packed cell volume which may be calculated as hemoglobin concentration (g/dl) / 34; this correction is needed because lithium is distributed in the plasma and not into the red or white cells on the first pass to the arterial circulation.

VALIDATION

Validation comparing the LiDCO System with bolus pulmonary artery catheter thermodilution technique (Figure 8) demonstrated a good overall agreement between the two methods (see Figure 8 $r^2 = 0.94$)¹¹. The conclusions were that a single bolus of lithium was at least as accurate as meaned triplicate bolus thermodilution. In another study where thermodilution and lithium dilution were compared to an aortic electromagnetic flow probe the LiDCO results showed less variability and therefore the LiDCO System was found to have a greater precision than single bolus thermodilution¹².

Larger animals (horses) and paediatric subjects have also been studied to ensure that the lithium dilution technique remained valid at extremes of flows (Figure 9). The body weight studied ranged from as small as a 2kg baby* up to a 550kg horse^{10,95}.




Figure 8. XY plot and Bland-Altman plot of 40 patients comparing LiDCO (average of 5 measurements) with bolus thermodilution (BTB) in 40 patients. XY plot: $\text{LiDCO} = 0.31 + 0.89 \text{ BTB}$ (L/minute) $r^2 = 0.94$. Bland-Altman plot: Mean differences BTB - LiDCO 0.25 L/minute, SD of the difference was 0.46 L/minute.




Figure 9. Comparisons ($n = 318$) of LiDCO vs bolus thermodilution in adults, paediatrics and horses.

*The use of LiDCOplus is unlicensed in patients <40kg (88lb).

LITHIUM CHLORIDE: THE FACTS

Multiple dosages of Lithium have been extensively investigated and the safety profile is well established. The pharmacokinetics of intravenous lithium chloride in man (and animals) has been documented⁷⁴, Lithium chloride has been used extensively in medicine for prophylactic and therapeutic treatment of unipolar and bipolar manic-depressive disorders^{75,76}. The lithium chloride is distributed throughout the total body water and excreted almost entirely by the kidneys. The half-life of lithium chloride in humans is 19.8 - 41.3 hrs^{77,78}. The recommended maximum total dose for a Lithium indicator dilution would have to be exceeded many times before toxic levels are reached. In fact, a single lithium chloride LiDCO indicator dilution determination at 0.3mmol is the equivalent to a steady state plasma lithium concentration of 1/240th of the therapeutic level. Lithium has been used for the measurement of cardiac output in thousands of patients over many years without a single side effect being reported.

(For warnings, contraindications and side-effect please refer to LiDCO Ltd Lithium Chloride Data Sheet and/or LiDCOplus manual)



- Erlanger & Hooker (1904) An experimental Study of blood pressure and of pulse pressure in man. *John Hopkins Hosp Rep* 12:145-378
- Pittman J, Bar Yosef S, SumPing J, Sherwood M, Mark J (2005) Continuous cardiac output monitoring with pulse contour analysis: A comparison with lithium indicator dilution cardiac output measurement. *Crit Care Med*. 33 (9) 2015-2021
- Aronson A, Heller L, Jayakar D, Jeevanandam V, Fisher M, Pfanzteller N, Babb P, Dupont F, Chaney M, Tung A. (2002) Continuous intraoperative cardiac output determination with arterial pulse wave analysis (PulseCO™) is valid and precise. Submitted to the Society of Cardiovascular Anesthesiologists Annual Meeting, May.
- Hamilton et al (2002) PulseCO: A less-invasive method to monitor cardiac output from arterial pressure after cardiac surgery. *Ann Thorac Surg* 74: 51408-12.
- Heller L B, Fisher M, Pfanzteller N, Jayakar D, Jeevanandam V, Aronson S. (2002) Continuous Intraoperative Cardiac Output Determination with Arterial Pulse Wave Analysis (PulseCO™) is Valid and Precise. *Anesth. Analg.* 2002;93:SCA1-SCA112
- Jonas M, Bruce R, Knight J, Kelly F, O'Brien T, Band D. (2001) Comparison of cardiac output measurements using a continuous arterial waveform analysis monitor (PulseCO) with an indicator dilution technique (LiDCO) in patients in intensive care. Accepted for presentation at the ICGM Sydney, October.
- Mason DJ, O'Grady M, Woods JP, McDonnell W. (2001) Assessment of lithium dilution cardiac output as a technique for measurement of cardiac output in dogs. *American Journal of Veterinary Research* Aug; 62(8): 1255 - 1261
- Mason DJ, O'Grady M, Woods JP, McDonnell W. (2000) Assessment of lithium chloride dilution as a new measure of cardiac output in the dog as compared to thermolulution. *Meeting of the American College of Veterinary Internal Medicine* May 2000.
- Linton RA, Jonas MM, Tibby SM, Murdoch IA, O'Brien TK, Linton NW, Band DM. (2000) Cardiac output measured by lithium dilution and transpulmonary thermolulution in patients in a paediatric intensive care unit. *Intensive Care Med*. Oct; 26(10): 1507-11.
- Linton RA, Young LE, Marlin DJ, Blissett KJ, Brearley JC, Jonas MM, O'Brien TK, Linton NW, Band DM Jones RS. (2000) Cardiac output measured by lithium dilution, thermolulution and transesophageal Doppler echocardiography in anesthetized horses. *American Journal of Veterinary Research* Jul; 61(7): 731-7.
- Linton R, Band D, O'Brien T, Jonas MM & Leach R. (1997) Lithium dilution cardiac output measurement: A comparison with thermolulution. *Crit Care Medicine*; 25: 1796-1800.
- Kuriita T, Morita K, Kato S, Kikura M, Horie M, Ikeda K. (1997) Comparison of the accuracy of the lithium dilution technique with the thermolulution technique for measurement of cardiac output. *British Journal of Anaesthesia*; 79: 770-775.
- Cohn JN, Pinkerson AL, Tristani FE: Mechanism of pulus paradoxus in clinical shock. *J. Clin. Invest.* 1967; 46:1744-55
- Reuter et al (2002) Optimizing fluid therapy in mechanically ventilated patients after cardiac surgery by on-line monitoring of left ventricular stroke variations. Comparison with aortic systolic pressure variations. *Brit Jour Anesth* 88 (1) 124-126
- Avila et al (2002) Predicting hypovolemia during mechanical ventilation: A prospective, clinical trial of doppler variations of aorta and axillary arterial velocities to identify systolic pressure variation. *Poster Presentation. American Society of Critical Care Anesthesiologists, Anesthesiology* 97 (3) 817
- Gunn and Pinsky (2001) Implications of arterial pressure variation in patients in the intensive care unit. *Critical Care* 7: 212-217
- Michard et al (1999) Clinical use of respiratory changes in arterial pulse pressure to monitor the hemodynamic effects of PEEP. *Critical Care Medicine* 159: 935-939
- Michard et al (2000) Relation between respiratory changes in arterial pulse pressure and fluid responsiveness in septic patients with acute circulatory failure. *Critical Care Medicine* 162: 134-138
- Michard & Teboul (2000) Using heart-lung interactions to assess fluid responsiveness during mechanical ventilation. *Crit Care* 4: 282-289
- Berkenstadt et al (2001) Stroke volume variation as a predictor of fluid responsiveness in patients undergoing brain surgery. *Anesth Analg* 92: 984-989
- Perel et al (1987) Systolic blood pressure variation is a sensitive indicator of hypovolemia in ventilated dogs subjected to graded hemorrhage. *Anesth* 67: 498-502
- Harrigan & Pinsky (2001) Heart-lung interactions. Part 1: effects of lung volume and ventilation as exercise. *Int Jour Inten Care Spring* 2001: 6-13
- Harrigan & Pinsky (2001) Heart-lung interactions. Part 2: effects of intrathoracic pressure. *Int Jour Inten Care Summer* 2001: 99-108
- Rooke et al (1995) The effect of graded hemorrhage and intravascular volume replacement on systolic pressure variation in humans during mechanical and spontaneous ventilation. *Anesth Analg* 80: 925-932
- Marik (1993) The systolic blood pressure variation as an indicator of pulmonary capillary wedge pressure in ventilated patients. *Anaes & Intens Care* 21 (4) 405-408
- Tavernier et al (1998) Systolic pressure variation as a guide to fluid therapy in patients with sepsis-induced hypotension. *Anesth* 89 (6) 1309-1310 [Abstract]
- Pizov et al (1996) Positive end-expiratory pressure-induced hemodynamic changes are reflected in the arterial pressure waveform. *Crit Care Med* 24 (8) 1381-7 [Abstract]
- Szold et al (1989) The effect of tidal volume and intravascular volume state on systolic pressure variation in ventilated dogs. *Inten Care Med* 15 (6) 368-371 [Abstract]
- Baeussier et al (1995) Determinants of systolic pressure variation in patients ventilated after vascular surgery. *J Cardioth Vasc Anesth* 9 (5) 547-551 [Abstract]
- Pizov et al (1990) The use of systolic pressure variation in hemodynamic monitoring during deliberate hypotension in spine surgery. *J Clin Anesth* 2 (2) 96-100 [Abstract]
- Perel et al (1987) Systolic blood pressure variation is a sensitive indicator of hypovolemia in ventilated dogs subjected to graded hemorrhage. *Anesth* 67 (4) 498-502 [Abstract]
- Weiss et al (1999) Systolic pressure variation in hemodynamic monitoring after severe blast injury. *J Clin Anesth* 11 (2) 132-135 [Abstract]
- Ornstein et al (1998) Systolic pressure variation predicts the response to acute blood loss. *J Clin Anesth* 10 (2) 137-140 [Abstract]
- Klingzang et al.(2002) Stroke volume variation as a predictor of fluid responsiveness for cardiac output in patients undergoing cardiac surgery. *SCCM* 173/M55 [Abstract]
- Reuter et al (2002) Left Ventricular stroke volume variations for functional preload monitoring after cardiac surgery in high risk patients. *SCCM* 93/M1[Abstract]
- R Boyd, RE Tremblay, Spencer FC et al (1959) Estimation of Cardiac output soon after intracardiac surgery with cardio-pulmonary bypass. *Ann Surg* 150: 613
- Clowes GHA, Del Guercio LRM (1960) Circulatory response to trauma of surgical operations. *Metabolism* 9:67-81
- Clowes et al, (1966) Early Study in Peritonitis. *Ann Surg* 1966:163; 866
- Shoemaker WC, Czer LS. (1979) Evaluation of the biologic importance of various Hemodynamic and oxygen transport variables: which should be monitored in postoperative shock? *Crit Care Med* 7: 424
- Shoemaker WC, Montgomery ES, Kaplan E, et al. (1973) Physiological patterns in surviving and non-surviving shock patients. Use of sequential cardiac respiratory variables in defining criteria for therapeutic goals and early warning of death. *Arch Surg* 106: 630
- Shoemaker WC, Pierchala C, Chang P et al (1977) Prediction of outcome and severity of illness by analysis of the frequency distributions of cardio respiratory variables. *Crit Care Med* 5: 82
- Shoemaker WC, Appel PL, Waxman K et al (1982) Clinical trial of survivors' cardio-respiratory patterns as therapeutic goals in critically ill postoperative patients. *Crit Care Med* 10: 398
- Edwards JD, Brown GC, Nightingale P, et al (1989) Use of survivors' cardio-respiratory values as therapeutic goals in septic shock. *Crit Care Med* 17:1098
- Schultz RJ, Whitfield GF, LaMura et al (1985) The role of physiologic monitoring in patients with fractures of the hip. *J Trauma* 25: 309
- Shoemaker et al (1988) Prospective trial of supranormal values of survivors as therapeutic goals in high-risk patients. *Chest* 94:1176-1186
- Berlauk et al (1991) Preoperative optimisation of cardiovascular hemodynamics improves outcome in peripheral vascular surgery: A prospective, randomised clinical trial. *Ann Surg* 214:289-297
- Fleming A, Bishop M, Shoemaker W et al (1992) Prospective trial of supranormal values as goals of resuscitation in severe trauma. *Arch Surg* 127:1175-1179
- Bishop MH, Shoemaker WC, Appel PL et al (1995) Prospective, randomised trial of survivor values of cardiac index, oxygen delivery, and oxygen consumption as resuscitation endpoints in severe trauma. *J Trauma*. 38:780-787
- Boyd O, Grounds RM, Bennett ED. (1993) A randomised clinical trial of the effect of deliberate perioperative increase of oxygen delivery on mortality in high-risk surgical patients. *JAMA* 270: 2699
- Wilson et al (1999) Reducing the risk of major elective surgery: randomised controlled trial of preoperative optimisation of oxygen delivery. *BMJ* 318:1099-1103.
- Sinclair et al (1997) Intraoperative intravascular volume optimisation and length of hospital stay after repair of proximal femoral fracture: randomised controlled trial. *BMJ* 315:909-912
- Shoemaker et al (1999) Intraoperative evaluation of tissue perfusion in high-risk patients by invasive and non-invasive Hemodynamic monitoring. *Crit Care Med* 27 (10) 2147-2152
- Shoemaker et al (1990) The efficacy of central venous and pulmonary artery catheters and therapy based upon them in reducing mortality and morbidity. *Arch Surg* 125:1332-1338
- Bland, Shoemaker et al (1985) Hemodynamic and oxygen transport patterns in surviving and non-surviving postoperative patients. *Crit Care Med* 13 (2) 85-90
- Gutierrez G, Palizas F, Doglio G et al (1992) Gastric intramucosal pH as a therapeutic index of tissue oxygenation in critically ill patients. *Lancet* 339: 19
- Grounds R.M., (2003) Reducing mortality and complications in patients undergoing surgery at high risk for post operative complications and death. In: *Recent Advances in Anaesthesia and intensive care* 22 ed: Adams, AP, Cashman, JN, Grounds, RM, Greenwich Medical Media Ltd, London.
- Grocott & Gan, Hemodynamic "Optimisation" Goal Is Improved Outcome
- Kern & Shoemaker (2002) Meta-analysis of hemodynamic optimisation in high-risk patients. *Crit Care Med* 30 (8) 1686-1692
- Bennett (2002) Goal-directed therapy is successful in the right patients. *Crit Care Med* 30 (8) 1909-10
- Singh and Manji (2001) A survey of pre-operative optimisation of high-risk surgical patients undergoing major elective surgery. *Anaesthesia* 56:988-1002
- Wilson et al (2001) Cardiovascular optimisation for high-risk surgery. *Cur Op Crit Care* 7:195-199
- Boyd & Bennett (1999) Achieving the goal. *Crit Care Med* 27:2298-9
- Nethrasisgamani & Fielden (2002) "Pre-operative optimisation of high risk surgical patients [Letter]" *Anaesthesia* 57:405-6
- Poeze et al (2000) Goal-orientated haemodynamic therapy: a plea for a closer look at using peri-operative oxygen transport. *Int Care Med* 26:635-637
- Traber (2002) Fluid resuscitation after hypovolemia. *Crit Care Med* 30 (8) 1922
- Parry-Jones AJD and Pittman JAL (2003) Arterial pressure and stroke volume variability as measurements for cardiovascular optimisation. *International Journal of Intensive Care, Summer* vol. 10 (2)
- Boyd, Grounds, Bennett (1993) The use of dopexamine hydrochloride to increase oxygen delivery perioperatively. *Anesth. Analg.* 76: 372-376
- Jonas M et al (2002) Real-time, continuous monitoring of cardiac output and oxygen delivery. *Inter J. Intensive Care* 9 (1) 33-42.
- Yu et al (1993) Affect of maximum oxygen delivery on morbidity and mortality in critically ill patients: a prospective, randomised controlled study. *Crit. Care Med* 21:830-838
- Rivers et al (2001) Early goal-directed therapy in the treatment of severe sepsis and septic shock. *N Engl J Med* 345 (19) 1368-77
- Hamilton et al (1932) Studies on the circulation IV. Further analysis of the injection methods and of changes in hemodynamics under physiological conditions. *Am.J. Physiol* 99:534-551
- Linton et al (1993) A new method of measuring cardiac output in man using lithium dilution. *Br.J. Anaesth.* 71:262-266
- Linton et al (1998) Lithium dilution cardiac output measurement – a brief review. In: Ikeda et al, *State of the art technology in anaesthesia & intensive care*. Elsevier 61-66.
- Jonas et al (2001) The pharmacokinetics of intravenous lithium chloride in patients and normal volunteers. *J. trace Elements Microprobe Technologies* 19:313-320
- Amidsen (1980) Lithium. In: Schentag et al, *Applied pharmacokinetics: principles of therapeutic drug monitoring*. San Francisco: Therapeutics Inc 586-617.
- Price et al (1994) Lithium in the treatment of mood disorders. *N Engl J Med*. 331:591-598
- Neilsen-Kudsk, Amidsen (1979) Analysis of the Pharmacokinetics in man. *Eur J Clin Pharmacol.* 16:271-277
- Mason, McQueen, Keary (1978) Pharmacokinetics of lithium elimination and half-life, renal clearance and apparent volume distribution in schizophrenia. *Clin Pharmacokinet* 3:241-246.
- Mappes et al (2001) Perioperative hemodynamic beat-to-beat monitoring of cardiac output:the new PulseCO system. *American Soc. Of Anesthesiology, New Orleans, Oct., Abstracts*
- Hamilton TT, Jessen ME., (2002) PulseCO: Accurate monitoring of cardiac output from the arterial waveform. Presented at the 8th Annual Current Trends in Thoracic Surgery VIII Meeting 2002.
- Teboul J-L (2003) Editorial: Dynamic concepts of volume responsiveness. *International Journal of Intensive Care, Summer* vol. 10 (2)
- Reuter et al (2002) Stroke volume variations for assessment of cardiac responsiveness to volume variations for assessment of cardiac responsiveness to volume loading in mechanically ventilated patients after cardiac surgery. *Intensive Care Med*. 28:392-398.
- Friesse, R., Lemaire, S., Minei, J. (2003) Reduction of pulse pressure variation is correlated with preload recruitable increases in cardiac index. University of Texas, Southwestern medical center, Dallas TX. *Poster abstract presented at the AAST (Tampa)*.
- Ziegler et al (1997) A prospective randomised trial of preoperative 'optimisation' of cardiac function in patients undergoing elective peripheral vascular surgery. *Surgery* 122:584-592.
- Polonen et al (2000) A prospective, randomised study of goal-oriented hemodynamic therapy in cardiac surgical patients. *Anaesth Analg* 90:1052-1059.
- Lobo et al (2000) Effects of maximizing oxygen delivery on morbidity and mortality in high-risk surgical patients. *Crit Care Med* 28:3396-3404.
- Smith J, Kirwan C, Lei K, Beale R (2005) Cardiac Output measured by lithium dilution and transpulmonary thermolulution in adult intensive care patients. *Crit Care Med*. 33(12) Suppl; 207-S; A56.
- Rupert Pearse, Deborah Dawson, Jayne Fawcett, Andrew Rhodes, R Michael Grounds and E David Bennett (2005) Early goal-directed therapy after major surgery reduces complications and duration of hospital stay. A randomised, controlled trial. [ISRCTN38797445], *Critical Care* 2005, 9:R687-R693 (DOI 10.1186/cc3887).
- Hallowell G, Carley K (2005) Use of lithium dilution and pulse contour analysis cardiac output determination in anesthetized horses: a clinical evaluation. *Vet Anaes & Analgesia* 32, 201-211.
- Costa MG, Lugano M, Coccia C, Della Rocca G (2004) Continuous and intermittent lithium dilution cardiac output measurements. *European J of Anaesthesiology* Vol. 21, Supplement 32, A-88.
- Costa MG, Cecconi M, Coccia C, Chiarandini P, Della Rocca G (2005) Continuous and intermittent cardiac output measurement, Pulmonary Artery Catheter vs. Lithium Dilution technique: Agreement at different levels of cardiac output. Presented at the Society of Critical Care Medicine, Phoenix Arizona, 15-19th January 2005
- Chinard FP, 'The Pulmonary Circulation and Interstitial Space' 1969 Ch6 79-98.
- Stewart GN, 'The Pulmonary Circulation time, the quantity of blood on the lungs and the output of the heart. *Am J Physiol* 1921. 58 20-44.
- Kirwan C et al (2005) A comparison of two calibrated, continuous, arterial pressure waveform based measurements of cardiac output over 24 hours. *Crit Care Med*. 33(12) Suppl; 207-S; A56.
- Kim JJ, Dreyer WJ, Chang AC, Breinholt JP, Grifka RG (2006) Arterial pulse wave analysis: An accurate means of determining cardiac output in children; *Pediatr Crit Care Med* 2006 Vol. 7, No. 6, 532-35.

Systolic Pressure Variation

Pulse Pressure Variation

Lithium Calibration

Absolute Cardiac Output

Global Sales:

Email: info@lidco.com

Telephone: +44 (0) 1223 830 666

Fax: +44 (0) 1223 837 241

LiDCO Ltd, Unit M, South Cambridge Business Park
Babraham Road, Sawston, Cambridge, CB22 3JH, UK


For USA Ordering:

Email: ussales@lidco.com

Tel: +1 877 LiDCO11 (877 543 2611)

Fax: +1 847 265 3737

LiDCO Ltd, 500 Park Ave., Suite103
Lake Villa, IL 60046, USA



Stroke Volume Variation