

Parameter	Equation	Normal Range
Normal Hemodynamic Parameters - Adult		
Arterial Blood Pressure (BP)	Systolic (SBP) Diastolic (DBP)	90 - 140 mmHg 60 - 90 mmHg
Mean Arterial Pressure (MAP)	$SBP + (2 \times DBP)/3$	70 - 105 mmHg
Systolic Pressure Variation (SPV)	$(SP_{max} - SP_{min})^\nabla$	<5 mmHg unlikely to be preload responsive >5mmHg likely to be preload responsive
Pulse Pressure Variation (PPV)%	$(PP_{max} - PP_{min}) / [(PP_{max} + PP_{min})/2] \times 100^\nabla$	<10% unlikely to be preload responsive >13-15% likely to be preload responsive
Stroke Volume Variation (SVV)%	$(SV_{max} - SV_{min}) / [(SV_{max} + SV_{min})/2] \times 100^\nabla$	<10% unlikely to be preload responsive >13-15% likely to be preload responsive
∇ = averaged over 10 sec. of BP data updated every 4 beats		
Right Atrial Pressure (RAP)		2 - 6 mmHg
Right Ventricular Pressure (RVP)	Systolic (RVSP) Diastolic (RVDP)	15 - 25 mmHg 0 - 8 mmHg
Pulmonary Artery Pressure (PAP)	Systolic (PASP) Diastolic (PADP)	15 - 25 mmHg 8 - 15 mmHg
Mean Pulmonary Artery Pressure (MPAP)	$[PASP + (2 \times PADP)]/3$	10 - 20 mmHg
Pulmonary Artery Wedge Pressure (PAWP)		6 - 12 mmHg
Left Atrial Pressure (LAP)		6 - 12 mmHg
Cardiac Output (CO)	$HR \times SV/1000$	4.0 - 8.0 l/min
Cardiac Index (CI)	CO/BSA	2.5 - 4.0 l/min/m ²
Stroke Volume (SV)	$CO/HR \times 1000$	60 - 100 ml/beat
Stroke Volume Index (SVI)	$CI/HR \times 1000$	33 - 47 ml/m ² /beat
Systemic Vascular Resistance (SVR)	$80 \times (MAP - RAP)/CO$	800 - 1200 dynes • sec/cm ⁵
Systemic Vascular Resistance Index (SVRI)	$80 \times MAP - RAP/CI$	1970 - 2390 dynes • sec/cm ⁵ /m ²
Pulmonary Vascular Resistance (PVR)	$80 \times (MPAP - PAWP)/CO$	<250 dynes • sec/cm ⁵
Pulmonary Vascular Resistance Index (PVRI)	$80 \times (MPAP - PAWP)/CI$	255 - 285 dynes • sec/cm ⁵ /m ²
Hemodynamic Parameters - Adult		
Left Ventricular Stroke Work (LVSW)	$SV \times (MAP - PAWP) \times 0.0136$	58 - 104 gm-m/beat
Left Ventricular Stroke Work Index (LVSWI)	$SVI \times (MAP - PAWP) \times 0.0136$	50 - 62 gm-m/m ² /beat
Right Ventricular Stroke Work (RVSW)	$SV \times (MPAP - RAP) \times 0.0136$	8 - 16 gm-m/beat
Right Ventricular Stroke Work Index (RVSWI)	$SVI \times (MPAP - RAP) \times 0.0136$	5 - 10 gm-m/m ² /beat
Coronary Artery Perfusion Pressure (CPP)	Diastolic BP - PAWP	60 - 80 mmHg
Right Ventricular End-Diastolic Volume (RVEDV)	SV/EF	100 - 160 ml
Right Ventricular End-Systolic Volume (RVESV)	EDV - SV	50 - 100 ml
Right Ventricular Ejection Fraction (RVEF)	SV/EDV	40 - 60%
Oxygenation Parameters - Adult		
Partial Pressure of Arterial Oxygen (PaO ₂)		80 - 100 mmHg
Partial Pressure of Arterial CO ₂ (PaCO ₂)		35 - 45 mmHg
Bicarbonate (HCO ₃)		22 - 28 mEq/l
pH		7.38 - 7.42
Arterial Oxygen Saturation (SaO ₂)		95 - 100%
Mixed Venous Saturation (SvO ₂)		60 - 80%
Arterial Oxygen Content (CaO ₂)	$(0.0138 \times Hgb \times SaO_2) + (0.0031 \times PaO_2)$	17 - 20 ml/dl
Venous Oxygen Content (CvO ₂)	$(0.0138 \times Hgb \times SvO_2) + (0.0031 \times PvO_2)$	12 - 15 ml/dl
A-V Oxygen Content Difference (C(a-v)O ₂)	$CaO_2 - CvO_2$	4 - 6 ml/dl
Oxygen Delivery (DO ₂)	$CaO_2 \times CO \times 10$	950 - 1150 ml/min
Oxygen Delivery Index (DO ₂ I)	$CaO_2 \times CI \times 10$	500 - 600 ml/min/m ²
Oxygen Consumption (VO ₂)	$(C(a - v)O_2) \times CO \times 10$	200 - 250 ml/min
Oxygen Consumption Index (VO ₂ I)	$(C(a - v)O_2) \times CI \times 10$	120 - 160 ml/min/m ²
Oxygen Extraction Ratio (O ₂ ER)	$[(CaO_2 - CvO_2)/CaO_2] \times 100$	22 - 30%
Oxygen Extraction Index (O ₂ EI)	$[SaO_2 - SvO_2]/SaO_2 \times 100$	20 - 25%